

Product datasheet for **HP216540**

FCRL1 Human qPCR Primer Pair (NM_052938)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	115350
Forward Sequence:	AGGCACAGACAATGGCGTCCAA
Reverse Sequence:	CAGATGAGGACCAGCCTGTCTC
ACCN:	NM_052938 , NM_052938.1 , NM_052938.2 , NM_052938.3 , NM_052938.4 , BC033690 , BC033690.1 , BC021274 , NM_052938.5
UniProt ID:	Q96LA6
Synonyms:	CD307a; FCRH1; IFGP1; IRTA5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)