

Product datasheet for **HP216583**

UGT2B28 Human qPCR Primer Pair (NM_053039)

Product data:

| | |
|-------------------|---|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 54490 |
| Forward Sequence: | GCCTTAGGTCTCAATACTCGGC |
| Reverse Sequence: | CTCATAGATGCCATTGGCTCCAC |
| ACCN: | NM_053039 , NM_053039.1 , BC156172 , BC157011 , NM_053039.2 |
| UniProt ID: | Q9BY64 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec. |
| Storage: | Store at -20 °C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)