

Product datasheet for **HP216605**

Dermcidin (DCD) Human qPCR Primer Pair (NM_053283)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	117159
Forward Sequence:	GGTTAGCCAGACAGGCACCAAA
Reverse Sequence:	CACGCTTCTAGATCTTCGACTG
ACCN:	<u>NM_053283</u> , <u>NM_053283.1</u> , <u>NM_053283.2</u> , <u>NM_053283.3</u> , <u>BC062682</u> , <u>BC062682.1</u> , <u>BC069108</u> , <u>BP313072</u> , <u>NM_053283.4</u>
UniProt ID:	<u>P81605</u>
Synonyms:	AIDD; DCD-1; DSEP; HCAP; PIF
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)