

Product datasheet for **HP216632**

PARD3B Human qPCR Primer Pair (NM_057177)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	117583
Forward Sequence:	CCAGTTGGTGATTCTGAGGAGG
Reverse Sequence:	TTCGCCTAGAGGTGGTTACAGC
ACCN:	<u>NM_057177</u> , <u>NM_057177.1</u> , <u>NM_057177.2</u> , <u>NM_057177.3</u> , <u>NM_057177.4</u> , <u>NM_057177.5</u> , <u>NM_057177.6</u> , <u>BC012984</u> , <u>BC156252</u> , <u>BC172499</u> , <u>BF434006</u> , <u>NM_057177.7</u>
UniProt ID:	<u>Q8TEW8</u>
Synonyms:	ALS2CR19; PAR3B; PAR3beta; PAR3L
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)