

## Product datasheet for **HP216642**

### ANTXR2 Human qPCR Primer Pair (NM\_058172)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	118429
Forward Sequence:	TCGGAATGGCAGTGTCTCTGTC
Reverse Sequence:	CAGGATAGGTGCAGGACAAAGC
ACCN:	<a href="#">NM_058172</a> , <a href="#">NM_058172.1</a> , <a href="#">NM_058172.2</a> , <a href="#">NM_058172.3</a> , <a href="#">NM_058172.4</a> , <a href="#">NM_058172.5</a> , <a href="#">BC034001</a> , <a href="#">BC042596</a> , <a href="#">BC047345</a> , <a href="#">BC107876</a> , <a href="#">BC114374</a> , <a href="#">BG235895</a> , <a href="#">BM989747</a> , <a href="#">BQ432519</a> , <a href="#">NM_058172.6</a>
UniProt ID:	<a href="#">P58335</a>
Synonyms:	CMG-2; CMG2; HFS; ISH; JHF
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)