

## Product datasheet for **HP216647**

### EVA1C Human qPCR Primer Pair (NM\_058187)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	59271
Forward Sequence:	CAACTCTCTGGCAGCCTTTGCT
Reverse Sequence:	TCCTTGGCACAGGACTCTCTGA
ACCN:	<a href="#">NM_058187</a> , <a href="#">NM_058187.1</a> , <a href="#">NM_058187.2</a> , <a href="#">NM_058187.3</a> , <a href="#">NM_058187.4</a> , <a href="#">BC038710</a> , <a href="#">BM836911</a> , <a href="#">BU175692</a> , <a href="#">BX342316</a> , <a href="#">BX373211</a> , <a href="#">NM_058187.5</a>
UniProt ID:	<a href="#">P58658</a>
Synonyms:	B18; B19; C21orf63; C21orf64; FAM176C; PRED34; SUE21
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)