

Product datasheet for **HP216667**

CLEC4D Human qPCR Primer Pair (NM_080387)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	338339
Forward Sequence:	GTAAGAGAGGCACAGGAGTGCA
Reverse Sequence:	GGAAGGCTCTCCAGTCAATAGG
ACCN:	NM_080387 , NM_080387.4 , BC032313 , BC032313.1 , BI909298 , NM_080387.5
UniProt ID:	Q8WXI8
Synonyms:	CD368; CLEC-6; CLEC6; CLECSF8; Dectin-3; MCL; MPCL
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)