

Product datasheet for **HP216670**

TCEAL2 Human qPCR Primer Pair (NM_080390)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	140597
Forward Sequence:	GGCAAGAGAGTTGAGGAACCGT
Reverse Sequence:	TGATCCTCCCTGCATCTCTGAC
ACCN:	NM_080390 , NM_080390.1 , NM_080390.3 , BC104764 , BC104764.1 , BC093932 , BM699243 , NM_080390.4
UniProt ID:	Q9H3H9
Synonyms:	my048; MY0876G05; WEX1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)