

Product datasheet for **HP216687**

B3GALT6 Human qPCR Primer Pair (NM_080605)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	126792
Forward Sequence:	ACCAGTACCTGGTGACGCACAA
Reverse Sequence:	GACCAGTCGTACCGTAGGACA
ACCN:	<u>NM_080605</u> , <u>NM_080605.1</u> , <u>NM_080605.2</u> , <u>NM_080605.3</u> , <u>BC015549</u> , <u>BC041621</u> , <u>BC051719</u> , <u>BC160034</u> , <u>BE249959</u> , <u>BF002665</u> , <u>BG290330</u> , <u>NM_080605.4</u>
UniProt ID:	<u>Q96L58</u>
Synonyms:	ALGAZ; beta3GalT6; EDSP2; EDSSPD2; SEMDJL1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)