

Product datasheet for **HP216692**

WFDC3 Human qPCR Primer Pair (NM_080614)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	140686
Forward Sequence:	TCGGATCTGCCGAGACATTCTT
Reverse Sequence:	CTACACAGCTCTTGTTCAGCC
ACCN:	BC026014 , NM_080614 , NM_080614.1 , BC102013 , BC102014 , BC146426 , BC148860 , BX091378 , NM_080614.2
UniProt ID:	Q8IUB2
Synonyms:	dj447F3.3; WAP14
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)