

Product datasheet for **HP216752**

OXGR1 Human qPCR Primer Pair (NM_080818)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	27199
Forward Sequence:	CTCGATGTGCAGTTGTAGCCTG
Reverse Sequence:	GAGTTCATCCGAAGCTGGTGAGG
ACCN:	<u>NM_080818</u> , <u>NM_080818.1</u> , <u>NM_080818.2</u> , <u>NM_080818.3</u> , <u>NM_080818.4</u> , <u>BC103881</u> , <u>BC103882</u> , <u>BC103883</u>
UniProt ID:	<u>Q96P68</u>
Synonyms:	aKGR; GPR80; GPR99; P2RY15; P2Y15
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)