

Product datasheet for **HP216768**

STK35 Human qPCR Primer Pair (NM_080836)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	140901
Forward Sequence:	AGCGCATGAGTCACGGCAACAA
Reverse Sequence:	GTCTCCACCTTCACAGAACTCC
ACCN:	<u>NM_080836</u> , <u>NM_080836.1</u> , <u>NM_080836.2</u> , <u>NM_080836.3</u> , <u>BC017340</u> , <u>BC111573</u> , <u>BC140883</u> , <u>BC140884</u> , <u>NM_080836.4</u>
UniProt ID:	<u>Q8TDR2</u>
Synonyms:	CLIK1; STK35L1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)