

Product datasheet for **HP216822**

MOBKL2A (MOB3A) Human qPCR Primer Pair (NM_130807)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	126308
Forward Sequence:	CGCAGATCAACAACGAGGACCT
Reverse Sequence:	GCCGTCAAAGTGGTGGATGTAG
ACCN:	<u>NM_130807</u> , <u>NM_130807.1</u> , <u>NM_130807.2</u> , <u>BC015049</u> , <u>BC015049.1</u> , <u>BC009547</u>
UniProt ID:	<u>Q96BX8</u>
Synonyms:	MOB-LAK; MOB1C; MOBKL2A; moblak
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)