

## Product datasheet for **HP216836**

### LRPPRC Human qPCR Primer Pair (NM\_133259)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10128
Forward Sequence:	ATCCGACATGGTTACTGGTGCC
Reverse Sequence:	GTGTCAAGGACAGCAGATGAATC
ACCN:	<u><a href="#">NM_133259</a></u> , <u><a href="#">NM_133259.1</a></u> , <u><a href="#">NM_133259.2</a></u> , <u><a href="#">NM_133259.3</a></u> , <u><a href="#">BC010282</a></u> , <u><a href="#">BC026034</a></u> , <u><a href="#">BC038181</a></u> , <u><a href="#">BC050311</a></u> , <u><a href="#">BC130285</a></u> , <u><a href="#">BQ021667</a></u>
UniProt ID:	<u><a href="#">P42704</a></u>
Synonyms:	CLONE-23970; GP130; LRP130; LSFC; MC4DN5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)