

Product datasheet for **HP216852**

GRIN3A Human qPCR Primer Pair (NM_133445)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	116443
Forward Sequence:	ACACGGCAAACCTTGCTGCTGT
Reverse Sequence:	CTTCAGCACTGCTTTCTCGGAC
ACCN:	<u>NM_133445</u> , <u>NM_133445.1</u> , <u>NM_133445.2</u> , <u>BC029757</u> , <u>BC132866</u> , <u>BC132868</u> , <u>BC167432</u>
UniProt ID:	<u>Q8TCU5</u>
Synonyms:	GluN3A; NMDAR-L; NMDAR3A; NR3A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)