

Product datasheet for **HP216869**

Phostensin (PPP1R18) Human qPCR Primer Pair (NM_133471)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	170954
Forward Sequence:	TACCCAACTGCCGAGGAGATCT
Reverse Sequence:	GTATTGGTACGTGGTCTCCAGG
ACCN:	<u>NM_133471</u> , <u>NM_133471.1</u> , <u>NM_133471.2</u> , <u>NM_133471.3</u> , <u>BC066644</u> , <u>BC006176</u> , <u>BC023314</u> , <u>BC047721</u> , <u>BX647622</u> , <u>NM_133471.4</u>
UniProt ID:	<u>Q6NYC8</u>
Synonyms:	HKMT1098; KIAA1949
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)