

## Product datasheet for **HP216882**

### DNA helicase HEL308 (HELQ) Human qPCR Primer Pair (NM\_133636)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	113510
Forward Sequence:	CTACTCCACAGGAGTGCTCTGT
Reverse Sequence:	CACGACCAGCTCTGCCAATCAT
ACCN:	<a href="#">NM_133636</a> , <a href="#">NM_133636.1</a> , <a href="#">NM_133636.2</a> , <a href="#">NM_133636.3</a> , <a href="#">BC011863</a> , <a href="#">BC015428</a> , <a href="#">BC065814</a> , <a href="#">BC095473</a> , <a href="#">BC111502</a> , <a href="#">BC141524</a> , <a href="#">BM986846</a> , <a href="#">NM_133636.4</a>
UniProt ID:	<a href="#">Q8TDG4</a>
Synonyms:	HEL308
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)