

Product datasheet for **HP216893**

DTX3L Human qPCR Primer Pair (NM_138287)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	151636
Forward Sequence:	CCAGGTTATGAGTCCTTTGGCAC
Reverse Sequence:	TGCAGTTCGCTGTATTCCAGGG
ACCN:	NM_138287 , NM_138287.1 , NM_138287.2 , NM_138287.3 , BC042191 , BC042191.1 , BC016846 , BC040372 , BC060509 , BX648267 , BX648645
UniProt ID:	Q8TDB6
Synonyms:	BBAP; RNF143
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)