

Product datasheet for **HP216906**

NALP6 (NLRP6) Human qPCR Primer Pair (NM_138329)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	171389
Forward Sequence:	ACTGTGCCATCTGAGCAGCCTC
Reverse Sequence:	TCACTGAGCCTGTTGTGGAGGA
ACCN:	NM_138329 , NM_138329.1 , NM_138329.2
UniProt ID:	P59044
Synonyms:	AVR; CLR11.4; NALP6; NAVR; NAVR/AVR; PAN3; PYPAF5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)