

Product datasheet for **HP216914**

CLEC12A Human qPCR Primer Pair (NM_138337)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	160364
Forward Sequence:	GCTCAGAATGCCAGCCTGTTGA
Reverse Sequence:	CCACGAGTGGAAATCTTCTCAGG
ACCN:	NM_138337 , NM_138337.1 , NM_138337.2 , NM_138337.3 , NM_138337.4 , NM_138337.5 , BC126289 , BC027967 , BC063424 , BC126291 , BE242761 , BQ017053
UniProt ID:	Q5QGZ9
Synonyms:	CD371; CLL-1; CLL1; DCAL-2; MICL
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)