

## Product datasheet for **HP216929**

### SHF Human qPCR Primer Pair (NM\_138356)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	90525
Forward Sequence:	CCCAGAGAAGGTCCTGAAAATG
Reverse Sequence:	CATACAGAGGCAAGGGCTGAGT
ACCN:	<a href="#">NM_138356</a> , <a href="#">NM_138356.1</a> , <a href="#">NM_138356.2</a> , <a href="#">BC007586</a> , <a href="#">BC101601</a> , <a href="#">BC101603</a> , <a href="#">BC143601</a> , <a href="#">BC171730</a> , <a href="#">BM669861</a> , <a href="#">BX332806</a> , <a href="#">BX416415</a> , <a href="#">BX446025</a>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)