

## Product datasheet for **HP216961**

### SP140L Human qPCR Primer Pair (NM\_138402)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	93349
Forward Sequence:	GGAAAGTGACCAAGCATGTGGC
Reverse Sequence:	TCTCGTCCCATCTGCTCCAG
ACCN:	<a href="#">BC004921</a> , <a href="#">NM_138402</a> , <a href="#">NM_138402.1</a> , <a href="#">NM_138402.2</a> , <a href="#">NM_138402.3</a> , <a href="#">NM_138402.4</a> , <a href="#">BC004921.1</a> , <a href="#">BC105007</a> , <a href="#">BC143624</a> , <a href="#">BQ232706</a> , <a href="#">BQ690359</a> , <a href="#">BX281497</a> , <a href="#">NM_138402.6</a>
UniProt ID:	<a href="#">Q9H930</a>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)