

Product datasheet for **HP216975**

GOLM2 Human qPCR Primer Pair (NM_138423)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	113201
Forward Sequence:	AAGCAGGCTACCAAGGACAGAG
Reverse Sequence:	GTCTGCCTCATACTCACCTACG
ACCN:	NM_138423 , NM_138423.1 , NM_138423.2 , NM_138423.3 , BC012124 , BC041786 , BC063480 , BG327802 , BI467176 , NM_138423.4
UniProt ID:	Q6P4E1
Synonyms:	CASC4; H63
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)