

## Product datasheet for **HP216989**

### **C6orf150 (MB21D1) Human qPCR Primer Pair (NM\_138441)**

#### **Product data:**

|                          |  |
|--------------------------|--|
| <b>Product Type:</b>     | qPCR Primer Pairs  |
| <b>Gene ID:</b>          | 115004   |
| <b>Forward Sequence:</b> | AGGAAGCAACTACGACTAAAGCC  |
| <b>Reverse Sequence:</b> | CGATGTGAGAGAAGGATAGCCG   |
| <b>ACCN:</b>             | <u><a href="#">NM_138441</a></u> , <u><a href="#">NM_138441.1</a></u> , <u><a href="#">NM_138441.2</a></u> , <u><a href="#">BC113606</a></u> , <u><a href="#">BC012928</a></u> , <u><a href="#">BC108714</a></u> , <u><a href="#">BC113608</a></u> , <u><a href="#">BC143694</a></u> , <u><a href="#">NM_138441.3</a></u>                    |
| <b>UniProt ID:</b>       | <u><a href="#">Q8N884</a></u>  |
| <b>Synonyms:</b>         | C6orf150; h-cGAS; MB21D1   |
| <b>Components:</b>       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.  |
| <b>Quality Control:</b>  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| <b>Storage:</b>          | Store at -20°C.  |
| <b>Stability:</b>        | The primer mix is stable for one year from date of shipping.   |



[View online »](#)