

Product datasheet for **HP216990**

CCDC124 Human qPCR Primer Pair (NM_138442)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	115098
Forward Sequence:	GCTGATGCCAAGAAGCAGAAGG
Reverse Sequence:	CGTCTCCTTCTTACGTTCCAGC
ACCN:	<u>NM_138442</u> , <u>NM_138442.1</u> , <u>NM_138442.2</u> , <u>NM_138442.3</u> , <u>BC013949</u> , <u>BC013949.2</u> , <u>BG829184</u> , <u>BM998230</u> , <u>BP309483</u> , <u>BQ669035</u>
UniProt ID:	<u>Q96CT7</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)