

## Product datasheet for **HP217044**

### **B3GNT6 Human qPCR Primer Pair (NM\_138706)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	192134
<b>Forward Sequence:</b>	CTACATGGGCATGTGTCTGGAG
<b>Reverse Sequence:</b>	CTAGCAGCAACTCGCGGTACAT
<b>ACCN:</b>	<u><a href="#">NM_138706</a></u> , <u><a href="#">NM_138706.1</a></u> , <u><a href="#">NM_138706.2</a></u> , <u><a href="#">NM_138706.3</a></u> , <u><a href="#">NM_138706.4</a></u> , <u><a href="#">BC103908</a></u> , <u><a href="#">BC025357</a></u> , <u><a href="#">BC103909</a></u> , <u><a href="#">BC103910</a></u>
<b>UniProt ID:</b>	<u><a href="#">Q6ZMB0</a></u>
<b>Synonyms:</b>	3-Gn-T6; B3Gn-T6; beta-1; beta3Gn-T6; BGnT-6
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)