

Product datasheet for **HP217079**

CCDC148 Human qPCR Primer Pair (NM_138803)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	130940
Forward Sequence:	CAGACAGCATCACACTTTGCAGC
Reverse Sequence:	CTCTGTTGCTCCAGCCTAAGTC
ACCN:	<u>BC021739</u> , <u>NM_138803</u> , <u>NM_138803.1</u> , <u>NM_138803.2</u> , <u>NM_138803.3</u> , <u>BC107597</u> , <u>BC015395</u> , <u>BC031994</u> , <u>BC105928</u> , <u>BC107596</u> , <u>NM_138803.4</u>
UniProt ID:	<u>Q8NFR7</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)