

Product datasheet for **HP217082**

CD200R (CD200R1) Human qPCR Primer Pair (NM_138806)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	131450
Forward Sequence:	CGTCTCCCATTTGACTGGCAAC
Reverse Sequence:	CCAAATGAATCCCACGATGGTCA
ACCN:	NM_138806 , NM_138806.1 , NM_138806.2 , NM_138806.3 , BC069661 , BC069661.1 , BC093890 , BC069721 , BC069743 , BC143393 , NM_138806.4
UniProt ID:	Q8TD46
Synonyms:	CD200R; HCRTR2; MOX2R; OX2R
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)