

Product datasheet for **HP217092**

PABIR3 Human qPCR Primer Pair (NM_138819)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	159091
Forward Sequence:	CTGAGAAGAGTCAACAGTGCCC
Reverse Sequence:	GATGGAACCATGAAAGGGAGGAG
ACCN:	<u>BC063491</u> , <u>BC065225</u> , <u>NM_138819</u> , <u>NM_138819.1</u> , <u>NM_138819.2</u> , <u>NM_138819.3</u> , <u>BC063491.1</u> , <u>BC065225.1</u> , <u>BC017868</u> , <u>BC017868.1</u> , <u>BC016784</u> , <u>BG532414</u> , <u>NM_138819.4</u>
UniProt ID:	<u>Q6P4D5</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)