

Product datasheet for **HP217116**

DEFB127 Human qPCR Primer Pair (NM_139074)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	140850
Forward Sequence:	CAAGGACATTGCAGGAAAATCTGC
Reverse Sequence:	GCTTTGGACGAGGTGGCTTTGT
ACCN:	<u>NM_139074</u> , <u>NM_139074.1</u> , <u>NM_139074.2</u> , <u>NM_139074.3</u> , <u>BC069486</u> , <u>BC109395</u> , <u>BC110085</u> , <u>NM_139074.4</u>
UniProt ID:	<u>Q9H1M4</u>
Synonyms:	bA530N10.2; C20orf73; DEF-27; DEFB-27; DEFB27; hBD-27
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)