

## Product datasheet for **HP217118**

### CCDC98 (FAM175A) Human qPCR Primer Pair (NM\_139076)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84142
Forward Sequence:	AGCTGCTCTACTCATCGACTGG
Reverse Sequence:	CCAGTGGACATACAGGAACCTG
ACCN:	<a href="#">BC016905</a> , <a href="#">NM_139076</a> , <a href="#">NM_139076.1</a> , <a href="#">NM_139076.2</a> , <a href="#">BC016905.1</a> , <a href="#">BC039573</a> , <a href="#">BC039573.1</a> , <a href="#">BQ775596</a> , <a href="#">BQ775930</a>
UniProt ID:	<a href="#">Q6UWZ7</a>
Synonyms:	ABRA1; CCDC98; FAM175A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)