

## Product datasheet for **HP217138**

### SLC39A11 Human qPCR Primer Pair (NM\_139177)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	201266
Forward Sequence:	AGGATCGCACTGCTCATCTTGG
Reverse Sequence:	GGCACTCTCAAAGGTAGCAGATG
ACCN:	<u>NM_139177, NM_139177.1, NM_139177.2, NM_139177.3, BC035631, BC035631.1, BU686998</u>
UniProt ID:	<u>Q8N1S5</u>
Synonyms:	C17orf26; ZIP-11; ZIP11
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)