

Product datasheet for **HP217139**

ALKBH3 Human qPCR Primer Pair (NM_139178)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	221120
Forward Sequence:	CCACTGCTAAGAGCCATCTCCA
Reverse Sequence:	TCAATCACTCGTGGCTCAGGAG
ACCN:	<u>NM_139178</u> , <u>NM_139178.1</u> , <u>NM_139178.2</u> , <u>NM_139178.3</u> , <u>BC015155</u> , <u>BC067257</u> , <u>BC103812</u> , <u>BC103813</u> , <u>BC103814</u> , <u>NM_139178.4</u>
UniProt ID:	<u>Q96Q83</u>
Synonyms:	ABH3; DEPC-1; DEPC1; hABH3; PCA1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)