

Product datasheet for **HP217145**

MTFMT Human qPCR Primer Pair (NM_139242)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	123263
Forward Sequence:	GCCGACTTTTGAATGAGGCTCTT
Reverse Sequence:	CCAGTAACTGTGTCTCCGTGAAG
ACCN:	<u>NM_139242</u> , <u>NM_139242.1</u> , <u>NM_139242.2</u> , <u>NM_139242.3</u> , <u>BC033687</u> , <u>BC016630</u> , <u>BQ774051</u>
UniProt ID:	<u>Q96DP5</u>
Synonyms:	COXPD15; FMT1; MC1DN27
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)