

Product datasheet for **HP217147**

Tomosyn (STXBP5) Human qPCR Primer Pair (NM_139244)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	134957
Forward Sequence:	GACTACCTCCAGAAAGCAGTGC
Reverse Sequence:	CCTCTGGAACAACAGTCCCTTC
ACCN:	<u>NM_139244</u> , <u>NM_139244.1</u> , <u>NM_139244.2</u> , <u>NM_139244.3</u> , <u>NM_139244.4</u> , <u>BC043173</u> , <u>BC043645</u> , <u>BC065713</u> , <u>BC107737</u> , <u>BC113382</u> , <u>BC113408</u>
UniProt ID:	<u>Q5T5C0</u>
Synonyms:	LGL3; LLGL3; Nbla04300
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)