

## Product datasheet for **HP217156**

### MCFD2 Human qPCR Primer Pair (NM\_139279)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	90411
Forward Sequence:	TCGCCACAAGAATTGCAGCTCC
Reverse Sequence:	CACTCATTAGTGGTGCCTGTTCA
ACCN:	<a href="#">NM_139279</a> , <a href="#">NM_139279.1</a> , <a href="#">NM_139279.2</a> , <a href="#">NM_139279.3</a> , <a href="#">NM_139279.4</a> , <a href="#">NM_139279.5</a> , <a href="#">BC040357</a> , <a href="#">BC040357.2</a> , <a href="#">BC037845</a> , <a href="#">BI553070</a> , <a href="#">NM_139279.6</a>
UniProt ID:	<a href="#">Q8NI22</a>
Synonyms:	F5F8D; F5F8D2; LMAN1IP; SDNSF
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)