

Product datasheet for **HP217157**

ORMDL3 Human qPCR Primer Pair (NM_139280)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	94103
Forward Sequence:	TAACCCACTGGGAGCAGATGGA
Reverse Sequence:	CTCATCAGGGACACGGTGTGA
ACCN:	<u>NM_139280</u> , <u>NM_139280.1</u> , <u>NM_139280.2</u> , <u>NM_139280.3</u> , <u>BC017087</u> , <u>BC017087.1</u> , <u>BC000638</u> , <u>BC018837</u> , <u>BC050710</u> , <u>BC054042</u> , <u>BC071833</u> , <u>BU633754</u> , <u>NM_139280.4</u>
UniProt ID:	<u>Q8N138</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)