

Product datasheet for **HP217180**

SFXN5 Human qPCR Primer Pair (NM_144579)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	94097
Forward Sequence:	TGTCATCAGCGCCGTCTCCATT
Reverse Sequence:	GAACGGCACAAACCTCTGGATG
ACCN:	<u>NM_144579</u> , <u>NM_144579.1</u> , <u>NM_144579.2</u> , <u>BC101311</u> , <u>BC013233</u> , <u>BC018640</u> , <u>BC046630</u> , <u>BC101312</u> , <u>BC101313</u> , <u>NM_144579.3</u>
UniProt ID:	<u>Q8TD22</u>
Synonyms:	BBG-TCC; SLC56A5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)