

Product datasheet for **HP217181**

C1orf85 (GLMP) Human qPCR Primer Pair (NM_144580)

Product data:

| | |
|--------------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 112770 |
| Forward Sequence: | CTGATGGTGCTCCCTAAGGACA |
| Reverse Sequence: | GGTCTTCCCAAAGGCTTTGCTG |
| ACCN: | NM_144580 , NM_144580.1 , NM_144580.2 , BC018757 , BC011575 , BC036340 , BM833126 , BP302818 , NM_144580.3 |
| UniProt ID: | Q8WWB7 |
| Synonyms: | C1orf85; NCU-G1 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)