

Product datasheet for **HP217194**

NIPA1 Human qPCR Primer Pair (NM_144599)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	123606
Forward Sequence:	GGTCCACGTTTCGTGCTACAGAA
Reverse Sequence:	CGCCGTGTAAGCCAGGAAGTTT
ACCN:	NM_144599 , NM_144599.1 , NM_144599.2 , NM_144599.3 , NM_144599.4 , BC025678 , BC103705 , BC156247 , BC172484 , BK001020 , BX648722
UniProt ID:	Q7RTP0
Synonyms:	FSP3; SLC57A1; SPG6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)