

## Product datasheet for **HP217225**

### FAM131A Human qPCR Primer Pair (NM\_144635)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	131408
Forward Sequence:	CGAACAAGAGGCTCGCTTTGCA
Reverse Sequence:	TCATCAGTGGAGTCCTCGCCAT
ACCN:	<u>NM_144635, NM_144635.1, NM_144635.2, NM_144635.3, NM_144635.4, BC037299, BC026221</u>
UniProt ID:	<u>Q6UXB0</u>
Synonyms:	C3orf40; FLAT715; PRO1378
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)