

## Product datasheet for **HP217233**

### IGJ (JCHAIN) Human qPCR Primer Pair (NM\_144646)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3512
Forward Sequence:	CCAGGATCATCCGTTCTTCCGA
Reverse Sequence:	CAAATCTGGTTCTCAATGGTGAGG
ACCN:	<u><a href="#">NM_144646</a></u> , <u><a href="#">NM_144646.1</a></u> , <u><a href="#">NM_144646.2</a></u> , <u><a href="#">NM_144646.3</a></u> , <u><a href="#">BC038982</a></u> , <u><a href="#">BC038982.1</a></u> , <u><a href="#">BP221304</a></u> , <u><a href="#">BP261870</a></u> , <u><a href="#">BP418109</a></u> , <u><a href="#">BX648421</a></u> , <u><a href="#">NM_144646.4</a></u>
UniProt ID:	<u><a href="#">P01591</a></u>
Synonyms:	IGCJ; IGJ; JCH
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)