

## Product datasheet for **HP217241**

### TCP10L Human qPCR Primer Pair (NM\_144659)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	140290
Forward Sequence:	AGCCAGGAAGAAATCTGCAGCG
Reverse Sequence:	GAGTGCCACTCTGATTCTTGTTGG
ACCN:	<a href="#">NM_144659</a> , <a href="#">NM_144659.1</a> , <a href="#">NM_144659.2</a> , <a href="#">NM_144659.3</a> , <a href="#">NM_144659.4</a> , <a href="#">NM_144659.5</a> , <a href="#">BC033855</a> , <a href="#">BC022024</a> , <a href="#">BC027970</a> , <a href="#">BG717680</a> , <a href="#">BG721020</a> , <a href="#">BI829783</a> , <a href="#">BX111775</a>
UniProt ID:	<a href="#">Q8TDR4</a>
Synonyms:	C21orf77; PRED77; TCP10A-1; TCP10A-2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)