

Product datasheet for **HP217287**

FAM82A1 (RMDN2) Human qPCR Primer Pair (NM_144713)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	151393
Forward Sequence:	AGTGAGTCTGGCAAGTCGGAGA
Reverse Sequence:	CCATAAGCACGAGCAAATCGCC
ACCN:	<u>BC024243</u> , <u>NM_144713</u> , <u>NM_144713.1</u> , <u>NM_144713.2</u> , <u>NM_144713.3</u> , <u>NM_144713.4</u> , <u>BQ017142</u> , <u>BR000689</u> , <u>BR000692</u> , <u>NM_144713.5</u>
UniProt ID:	<u>Q96LZ7</u>
Synonyms:	BLOCK18; FAM82A; FAM82A1; PRO34163; PYST9371; RMD-2; RMD2; RMD4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)