

Product datasheet for **HP217303**

SMCR8 Human qPCR Primer Pair (NM_144775)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	140775
Forward Sequence:	GCTAAGCCCGTGAAACATTGGG
Reverse Sequence:	TGTC AAGGTCCAGGATGCTCGT
ACCN:	NM_144775 , NM_144775.1 , NM_144775.2 , BC142680 , BC001018 , BC005067 , BC014179 , BC040699 , BC101115 , BC101116 , BC101117 , BP280386 , NM_144775.3
UniProt ID:	Q8TEV9
Synonyms:	DENND8A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)