

## Product datasheet for **HP217338**

### NLRP11 Human qPCR Primer Pair (NM\_145007)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	204801
Forward Sequence:	CTGGATTACTGCTGTACCTGAG
Reverse Sequence:	CTCTCCAGTAGACAAGGCTCTTC
ACCN:	<u><a href="#">NM_145007</a></u> , <u><a href="#">NM_145007.1</a></u> , <u><a href="#">NM_145007.2</a></u> , <u><a href="#">NM_145007.3</a></u> , <u><a href="#">BC034730</a></u> , <u><a href="#">BC034730.1</a></u> , <u><a href="#">BC036506</a></u> , <u><a href="#">BU622142</a></u>
UniProt ID:	<u><a href="#">P59045</a></u>
Synonyms:	CLR19.6; NALP11; NOD17; PAN10; PYPAF6; PYPAF7
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)