

Product datasheet for **HP217352**

SPATS1 Human qPCR Primer Pair (NM_145026)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	221409
Forward Sequence:	CCTCTGGCAAAAGTGTCAGTTCC
Reverse Sequence:	GAGAGTTCCGATCCAAGTCAGC
ACCN:	NM_145026 , NM_145026.2 , NM_145026.3 , BC100956 , BC100957 , BC100958 , BC100959 , BC110065 , BC111414 , NM_145026.4
UniProt ID:	Q496A3
Synonyms:	DDIP; SPATA8; SRSP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)