

## Product datasheet for **HP217366**

### ODAD3 Human qPCR Primer Pair (NM\_145045)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	115948
Forward Sequence:	CCGAGCACATTACCAGCGTGTA
Reverse Sequence:	CTCCAGCTCATGTTTGGTCTC
ACCN:	<u><a href="#">BC014252</a></u> , <u><a href="#">NM_145045</a></u> , <u><a href="#">NM_145045.1</a></u> , <u><a href="#">NM_145045.2</a></u> , <u><a href="#">NM_145045.3</a></u> , <u><a href="#">NM_145045.4</a></u> , <u><a href="#">BC142637</a></u> , <u><a href="#">BC036041</a></u> , <u><a href="#">BC050355</a></u> , <u><a href="#">BC141828</a></u> , <u><a href="#">BC144421</a></u> , <u><a href="#">BI548500</a></u> , <u><a href="#">NM_145045.5</a></u>
UniProt ID:	<u><a href="#">A5D8V7</a></u>
Synonyms:	CCDC151; CILD30; ODA10
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)