

Product datasheet for **HP217383**

TRPV3 Human qPCR Primer Pair (NM_145068)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	162514
Forward Sequence:	ATCCTACTGCGGAGTGGCAACT
Reverse Sequence:	CGCTTCTCCTTGATCTCACGAC
ACCN:	NM_145068 , NM_145068.1 , NM_145068.2 , NM_145068.3 , BC104866 , BC104866.1 , BC104868 , BC143294 , BC143295 , BC143299 , BC143330 , BC143333 , BC143334 , BC143336 , BX537539 , NM_145068.4
UniProt ID:	Q8NET8
Synonyms:	FNEPPK2; OLMS; OLMS1; VRL3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)