

Product datasheet for **HP217385**

FAM200A Human qPCR Primer Pair (NM_145111)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	221786
Forward Sequence:	CATCGCCACAAGTTCTCAGTCG
Reverse Sequence:	TTTTCAGCCGCTGTGTGAGCCA
ACCN:	<u>NM_145111</u> , <u>NM_145111.1</u> , <u>NM_145111.2</u> , <u>NM_145111.3</u> , <u>BC040214</u> , <u>BC040214.1</u> , <u>BC009536</u> , <u>BC126337</u> , <u>BC130371</u> , <u>NM_145111.4</u>
UniProt ID:	<u>Q8TCP9</u>
Synonyms:	C7orf38
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)